



0.3mm-0.7mm Feldspathic EMAX Ceramic Veneers Natural Aesthetic Emax Dental Veneers

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Application: Suitable For Anterior And Posterior Teeth
- Esthetics: Highly Aesthetic And Natural Appearance
- Surface Texture: Smooth
- Thickness: 0.3-0.7mm
- Customization: Customizable Shape And Size
- Material: Feldspathic Ceramic
- Translucency: High Translucency
- Durability: Long-lasting
- Highlight: **0.7mm Feldspathic EMAX Ceramic Veneers, 0.3mm EMAX Ceramic Veneers, Natural Aesthetic Emax Dental Veneers**



More Images



Product Description

High Translucency 0.3-0.7mm Feldspathic Ceramic Veneer with Natural Aesthetic Appearance

Product Description:

These veneers are made from high-quality feldspathic porcelain, which is known for its durability and natural tooth color. With a thickness of 0.3-0.7mm, they are thin enough to be minimally invasive, requiring only a small amount of preparation before application. One of the key benefits of our Feldspathic Veneers is their high translucency, which allows them to blend seamlessly with your natural teeth for a more natural-looking smile. And because they are made from porcelain, they are long-lasting and resistant to staining and wear.

If you're looking for a non-invasive solution to improve the appearance of your teeth, our Feldspathic Porcelain Veneers may be just what you need. Contact us today to learn more!

Features:

Product Name: Feldspathic Veneer

Thickness: 0.3-0.7mm

Translucency: High Translucency

Application: Suitable For Anterior And Posterior Teeth

Surface Texture: Smooth

Color: Natural Tooth Color

Feldspathic Veneer, Feldspathic Porcelain Veneer, and Feldspathic Porcelain Veneers are all terms used to describe this product.

Technical Parameters:

Esthetics	Highly Aesthetic And Natural Appearance
Material	Feldspathic Ceramic
Color	Natural Tooth Color
Surface Texture	Smooth
Translucency	High Translucency
Feature	Non-invasive And Minimal Prep Is Required
Durability	Long-lasting
Customization	Customizable Shape And Size
Application	Suitable For Anterior And Posterior Teeth
Thickness	0.3-0.7mm

Applications:

The highly aesthetic and natural appearance of Feldspathic Veneers makes them suitable for both anterior and posterior teeth, providing a seamless and natural-looking solution for a variety of dental problems. Their natural tooth color ensures that they blend in seamlessly with the patient's existing teeth, providing a uniform and consistent smile.

The thickness of these veneers ranges from 0.3-0.7mm, providing dentists with the flexibility to choose the appropriate thickness for their patients' needs. These veneers are packaged in plastic bags and box, ensuring safe and secure delivery to the dentist's office.

ADL offers a negotiable minimum order quantity, price, and delivery time, providing dentists with the flexibility to choose the right order for their needs. Payment can be made through Paypal or T/T, making the ordering process simple and straightforward.

With a supply of 630 technicians, ADL's Feldspathic Veneers are a reliable and accessible solution for dental professionals. These veneers can be used in a variety of occasions and scenarios, providing patients with a natural-looking and highly aesthetic solution to their dental problems.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China

