



Customized Titanium Base Corrosion Resistance Dental Implants And Solutions

Our Product Introduction

for more products please visit us on adldentalab.com


Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Customization: Can Be Customized According To Specific Requirements
- Corrosion Resistance: High Corrosion Resistance
- Material: Titanium
- Biocompatibility: Biocompatible With Human Body
- Precision: Precisely Machined For Accurate Fit
- Compatibility: Compatible With Most Dental Implants
- Size: Customizable
- Color: Silver And Gold
- Highlight: **Customized Titanium Base, Titanium Base Corrosion Resistance, Titanium Base Dental Implants And Solutions**



More Images



Product Description

Customized Titanium Base for Corrosion Resistance Dental Implant Compatibility

Product Description:

One of the key features of the Customized Titanium Base is its ability to be customized according to specific requirements. This means that it can be tailored to fit the needs of each individual patient, ensuring that it provides the best possible fit and comfort. Customization options include size, shape, and design, making it a versatile and highly adaptable product.

The Customized Titanium Base is also biocompatible with the human body, which means that it is safe to use and poses no risk of adverse reactions or side effects. This is an essential feature for any dental implant product, as it ensures that patients are not exposed to any harmful substances or materials.

When it comes to dental implants, the quality and durability of the base is of utmost importance. The Customized Titanium Base is made from high-grade titanium, which is known for its strength, durability, and corrosion resistance. This ensures that the base will last for many years and provides a stable foundation for the implant.

The Customized Titanium Base is an excellent choice for patients who are looking for a high-quality, customizable, and biocompatible dental implant base. Its compatibility with most dental implants, its customizable features, and its biocompatibility make it a popular choice among dentists and dental technicians. If you are looking for a reliable and long-lasting dental implant base, the Customized Titanium Base is the perfect product for you.

Features:

Product Name: Customized Titanium Base

Size: Customizable

Biocompatibility: Biocompatible With Human Body

Compatibility: Compatible With Most Dental Implants

Color: Silver And Gold

Precision: Precisely Machined For Accurate Fit

Technical Parameters:

Precision	Precisely Machined For Accurate Fit
Color	Silver And Gold
Corrosion Resistance	High Corrosion Resistance
Customization	Can Be Customized According To Specific Requirements
Biocompatibility	Biocompatible With Human Body
Size	Customizable
Material	Titanium
Compatibility	Compatible With Most Dental Implants

Applications:

The Customized Titanium Base is compatible with most dental implants, making it perfect for dental implant procedures. Its high corrosion resistance makes it durable and long-lasting. This product can be customized according to specific requirements, ensuring that it will fit perfectly with your dental implant system.

The Customized Ti-base is biocompatible with the human body, which means that it will not cause any harm or negative reaction to the patient's body. This is important because it ensures that the implant is safe and will not cause any complications or infections.

The Customized Titanium Base is suitable for a variety of dental implant applications, including single-tooth implants, multiple-tooth implants, and full-arch implants. It can also be used for both upper and lower jaw implants, making it a versatile product.

Whether you are a dentist or a dental technician, the Customized Titanium Base is an excellent choice for your dental implant needs. It is easy to use, durable, and customizable, making it a reliable and cost-effective option.

With a supply ability of 630 technicians, ADL's Customized Titanium Base is readily available and can be purchased using Paypal or T/T. So, if you are in need of a high-quality and customizable titanium base for your dental implant procedures, look no further than ADL.



Advanced Dental Laboratory

+86 18665847934

info@adldentalab.com

adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong

