



Internal Hex Custom Titanium Abutment Screw Cement Retained Varying Shape

Our Product Introduction

for more products please visit us on adldentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians




Product Specification

- Material: Titanium
- Shape: Varies Based On Customization
- Connection Type: Internal Hex
- Customization Options: Abutment Height, Diameter, Angulation, Shape, Margin Design, Retention Mechanism, Prosthetic Platform, Prosthetic Type
- Color: Silver And Gold
- Diameter: Varies Based On Customization
- Surface Treatment: Polished
- Retention Mechanism: Screw-retained Or Cement-retained
- Highlight: **Internal Hex Custom Titanium Abutment, Cement Retained Titanium Abutment, Screw Retained titanium abutments dental**



More Images



Product Description

Custom Titanium Abutment Internal Hex Screw Cement Retained Varying Shape / Color for Customization

Product Description:

The Customize Titanium Abutment offered by our company is a premium quality dental product that is designed to provide a customized solution for your patients. Our titanium custom abutment is created to provide maximum stability and durability, while also offering a high degree of customization to meet the unique needs of your patients.

Our titanium custom abutment offers a wide range of customization options, including abutment height, diameter, angulation, shape, margin design, retention mechanism, prosthetic platform, and prosthetic type. This allows you to create a tailored solution for each patient, ensuring that their needs are met and that they receive the best possible care.

The prosthetic platform of our titanium custom abutment can be narrow or wide, depending on the needs of your patient. This allows you to create a more stable and secure fit for their prosthetic, ensuring that it stays in place and does not shift or move around.

The angulation of our titanium custom abutment varies based on customization, allowing you to create a more natural and comfortable fit for your patient. This helps to reduce discomfort and ensure that your patient is able to speak, eat, and drink normally.

The retention mechanism of our titanium custom abutment can be either screw-retained or cement-retained, depending on your preferences and the needs of your patient. This allows you to create a secure and reliable fit for their prosthetic, ensuring that it stays in place and does not become dislodged or damaged.

The height of our titanium custom abutment varies based on customization, allowing you to create a more comfortable and natural fit for your patient. This helps to reduce discomfort and ensure that your patient is able to speak, eat, and drink normally.

Overall, our Customize Titanium Abutment is an excellent choice for dental professionals who want to provide the highest quality care to their patients. With a wide range of customization options and a focus on durability and stability, our titanium custom abutment is designed to meet the unique needs of each patient and ensure that they receive the best possible care. So, if you're looking for a reliable and customizable solution for your patients, be sure to check out our titanium custom abutment today!

Features:

Product Name: Customize Titanium Abutment

Height: Varies Based On Customization

Compatibility: Compatible With Most Implant Systems

Prosthetic Type: Single Or Multi-unit

Prosthetic Platform: Narrow Or Wide

Shape: Varies Based On Customization

Technical Parameters:

Technical Parameter	Value
Prosthetic Platform	Narrow or Wide
Height	Varies Based on Customization
Angulation	Varies Based on Customization
Compatibility	Compatible with Most Implant Systems
Material	Titanium
Prosthetic Type	Single or Multi-unit
Connection Type	Internal Hex
Margin Design	Chamfer or Shoulder
Customization Options	Abutment Height, Diameter, Angulation, Shape, Margin Design, Retention Mechanism, Prosthetic Platform, Prosthetic Type
Diameter	Varies Based on Customization

Applications:

The Customize Titanium Abutment is a custom implant Ti abutment that is used to support prosthetics such as a dental crown, bridge, or denture. This product can be used for single or multi-unit prosthetic type and is available in narrow or wide prosthetic platforms. The shape of the abutment varies based on customization, ensuring a perfect fit for the patient's needs.

The Customize Titanium Abutment is an ideal solution for patients who require a customized implant solution. The product is made of titanium, a biocompatible material that is widely used in the medical industry. The surface treatment of the abutment is polished, ensuring a smooth and comfortable fit for the patient. With a minimum order quantity that is negotiable and customizable, this product is perfect for dental clinics and laboratories of any size and budget.

The Customize Titanium Abutment is packaged in plastic bags and boxes, ensuring its safety during transportation. The delivery time for this product is only 72 hours, making it an ideal solution for patients who require immediate implant placement. Payment can be made using Paypal or T/T, ensuring a convenient and hassle-free transaction. With a supply ability of 630 technicians, ADL can meet the demand for the Customize Titanium Abutment.

In conclusion, the Customize Titanium Abutment is an excellent product that is suitable for a wide variety of dental patients. This product is customizable, ensuring a perfect fit for the patient's needs. With its biocompatible material, polished surface treatment, and quick delivery time, the Customize Titanium Abutment is an ideal solution for dental clinics and laboratories worldwide.



Advanced Dental Laboratory



+86 18665847934



info@adldentalab.com



adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China