



Silver Gold Titanium Multi Unit Abutment For Single Or Multi Unit Prosthetic Platform

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Material: Titanium
- Shape: Varies Based On Customization
- Color: Silver And Gold
- Prosthetic Platform: Narrow Or Wide
- Prosthetic Type: Single Or Multi-unit
- Compatibility: Compatible With Most Implant Systems
- Customization Options: Abutment Height, Diameter, Angulation, Shape, Margin Design, Retention Mechanism, Prosthetic Platform, Prosthetic Type
- Diameter: Varies Based On Customization
- Highlight: **Silver Gold Titanium Multi Unit Abutment, Silver Gold Titanium Multi Unit Abutment cho, Silver Gold Multi Unit Abutment cho**



More Images



Product Description

Silver and Gold Custom Titanium Multi-unit Abutment for Single or Multi-unit Prosthetic Platform

Product Description:

One of the main advantages of the Custom Titanium Multi-unit Abutment is its customization options. Abutment height, diameter, angulation, shape, margin design, retention mechanism, prosthetic platform, and prosthetic type can all be tailored to meet the unique needs of the patient.

The prosthetic platform of the Multi-unit can be either narrow or wide, depending on the patient's needs. This allows for greater flexibility in choosing the right prosthetic type for each patient.

Customizing the abutment height, diameter, and angulation is essential for creating a precise fit for the final restoration. The margin design can also be customized to ensure that the abutment meets the patient's aesthetic requirements.

The retention mechanism is an important factor in ensuring that the abutment stays securely in place. The Multi-unit offers customizable retention mechanisms, ensuring that each patient's needs are met.

Overall, the Multi-unit is an excellent choice for creating custom aesthetic abutments. Its high-quality titanium construction, customization options, and ability to meet the unique needs of each patient make it an essential component in providing high-quality dental care.

Features:

Product Name: Custom Titanium Multi Unit Abutment

Prosthetic Type: Single Or Multi-unit

Height: Varies Based On Customization

Connection Type: Internal Hex

Margin Design: Chamfer Or Shoulder

Angulation: Varies Based On Customization

Custom Titanium Abutments are the perfect solution for those who require a unique fit and function. With the Customize Titanium Abutment, you can choose from a variety of customization options to create a custom implant Ti abutment that meets your specific needs. Whether you need a single or multi-unit prosthetic type, varying heights, internal hex connection type, or chamfer or shoulder margin design, our Customize Titanium Abutment can accommodate your needs. Contact us today to learn more about how we can help you create a custom implant Ti abutment that fits perfectly.

Technical Parameters:

Product Name	Custom Titanium Multi-unit Abutment
Prosthetic Type	Single Or Multi-unit
Material	Titanium
Color	Silver And Gold
Connection Type	Internal Hex
Angulation	Varies Based On Customization
Margin Design	Chamfer Or Shoulder
Diameter	Varies Based On Customization
Height	Varies Based On Customization
Compatibility	Compatible With Most Implant Systems
Retention Mechanism	Screw-retained Or Cement-retained

Applications:

Custom Titanium Multi-unit Abutment provides a retention mechanism of screw-retained or cement-retained. The abutments are made of pure titanium material, which is biocompatible, non-corrosive, and has excellent mechanical properties. The angulation of the abutment varies based on customization, and customization options include abutment height, diameter, angulation, shape, margin design, retention mechanism, prosthetic platform, and prosthetic type.

The prosthetic platform of the Custom Titanium Multi-unit Abutment is available in narrow or wide options. The customization options allow for a perfect fit to the patient's anatomy, creating a natural and comfortable tooth replacement. The Custom Titanium Multi-unit Abutment are suitable for various occasions and scenarios, including single tooth replacements, multiple missing teeth, and full-arch restorations.

The screw-retained retention mechanism provides easy access for the dentist to remove the abutment and prosthetic for maintenance and cleaning. The cement-retained retention mechanism provides excellent esthetics, as the abutment and prosthetic are cemented in place and appear as a natural tooth.

In conclusion, ADL's Custom Titanium Multi-unit Abutment offer a high-quality and customizable dental implant solution for various occasions and scenarios. The abutments have a retention mechanism of screw-retained or cement-retained, and the material is pure titanium, which is biocompatible and non-corrosive. With customization options such as abutment height, diameter, angulation, shape,


margin design, retention mechanism, prosthetic platform, and prosthetic type, the Custom Titanium Multi-unit Abutment provides a perfect fit to the patient's anatomy, creating a natural and comfortable tooth replacement.



Advanced Dental Laboratory

 +86 18665847934

 info@adldentalab.com

 adldentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China