



# CE FDA Ball Attachment Dentures Comfortable Denture Ball Attachment Long Lasting

Our Product Introduction

for more products please visit us on [adidentalab.com](http://adidentalab.com)

## Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians




## Product Specification

- Compatibility: Compatible With Dental Implants
- Ease Of Use: Easy To Insert And Remove
- Stability: Provides Stability For Dentures
- Retention: Secure Retention Of Dentures
- Comfort: Comfortable For The Patient
- Durability: Long-lasting
- Material: Titanium
- Adjustability: Adjustable For Optimal Fit
- Highlight: **fda ball attachment dentures, comfortable denture ball attachment, ball attachment dentures comfortable**



## More Images



## Product Description

### Comfortable Ball Attachment for Dentures Secure Stability And Durability

#### Product Description:

The Ball Attachment product is easy to insert and remove, making it a convenient choice for patients who require denture retention. This product is compatible with dental implants, which means that it can be used with any implant-supported denture system.

The one-piece ball attachments are made of high-quality materials, ensuring that they are durable and long-lasting. The natural-looking aesthetics of this product make it an ideal choice for patients who want to enhance the appearance of their dentures.

The Ball Attachment product provides secure retention of dentures, which is essential for patients who want to avoid embarrassment or discomfort caused by loose dentures. The implant ball attachment overdenture system is a popular choice for patients who have dental implants and want to improve the stability and comfort of their dentures.

In summary, the Ball Attachment product is an essential component of the implant ball attachment overdenture system. It provides secure retention of dentures, is easy to insert and remove, and is compatible with dental implants. With various sizes available, patients can choose the right size for their needs. The natural-looking aesthetics of this product make it a popular choice for patients who want to enhance the appearance of their dentures.

#### Features:

**Product Name: Ball Attachment**

**Material: Titanium**

**Shape: Round**

**Adjustability: Adjustable For Optimal Fit**

**Retention: Secure Retention Of Dentures**

**Comfort: Comfortable For The Patient**

**Keywords: Implant Ball Attachment Overdenture, One-Piece Ball Attachments, Implant Ball Attachment Overdenture**

#### Technical Parameters:

Material:	Titanium
Comfort:	Comfortable For The Patient
Compatibility:	Compatible With Dental Implants
Aesthetics:	Natural-looking
Shape:	Round
Hygiene:	Easy To Clean And Maintain
Retention:	Secure Retention Of Dentures
Ease Of Use:	Easy To Insert And Remove
Durability:	Long-lasting
Size:	Various Sizes Available

#### Applications:

Our Ball Attachment Overdenture is also long-lasting, providing you with durable support for your dentures. The adjustable nature of the product ensures that you can achieve an optimal fit, making it comfortable to wear. The natural-looking aesthetics of the Ball Attachment Overdenture will also give you the confidence to smile without worrying about the appearance of your dentures.

Our Ball Attachment Overdenture is also easy to clean and maintain, ensuring that you can maintain good hygiene and prevent the buildup of bacteria. The product is delivered in plastic bags and boxes, and we offer a fast delivery time of 72 hours.

Our Ball Attachment Overdenture is perfect for a range of application occasions and scenarios, including:

**Patients who have lost all or most of their teeth and require implant-supported dentures**

**Patients who want to improve the retention and stability of their dentures**

**Patients who want to improve the aesthetics of their dentures**

**Dental clinics and practitioners who want to offer a high-quality, reliable product to their patients**

We offer a negotiable minimum order quantity and price, making it accessible for a range of patients and dental practitioners. Payment can be made through Paypal or T/T, and we have a supply ability of 630 technicians.

Invest in the Ball Attachment Overdenture from ADL today and experience the benefits of a reliable, durable, and natural-looking product that will improve the retention and aesthetics of your dentures.



**Advanced Dental Laboratory**



+86 18665847934



info@adidentalab.com



adidentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China